

# User Guide of CpG\_MPs v1.1.0

CpG\_MPs is developed by Perl script which could be run in Linux or Window system. Users need to install Perl in the Window system before running the program. CpG\_MPs mainly includes four modules including Data Normalization, Identification of Methylation Patterns of regions, Identification of Differentially and Conserved Methylated Regions, and Calculation the sequence features of the regions of DNA methylation patterns.

## (1) Data Normalization

This module help you to convert methylation files into ".wig" format to display in ucsc.".wig" format is also the normal format which CpG\_MPs requires.

-----  
USAGE: Normalization [OPTION] {--m\_u <c1,c2>|--r\_t <c1,c2>} [-i <file|folder>] [-o <folder>]  
-----

### ARGUMENTS:

-h,--help,--man

Display the helpfile.

### OPTION:

--depth

A depth threshold,DEFAULT:5.A site is valid when its depth is more than the threshold.

--head

The methy files has a header.

(The command requires only one in these following two options.)

-m\_u,--methy\_unmethy

The column of methylated counts and unmethylated counts.

-r\_t,--ratio\_total

The column of methylation level(from 0 to 1) and total counts(i.e. read depth).

### I/O:

-i,--infile

The input methy file.Also can be a folder which has methy files splited by chromosome.REQUIRED

-o,--outdir

The output folder.REQUIRED

-----  

### EXAMPLE:

#### BRIEF:

Normalization --methy\_unmethy 5,6 -i test.txt -o out\_folder

Split "test.cov" into "out\_folder/" by chromosomes. methy counts is in col 5, unmethy counts is in col 6.

(the format of test.txt seems like ".cov" files comes from Bismark .)

#### OUTPUT:

Wig format will like this (chr1.wig):

```
track type=wiggle_0 name="Methylation Level" visibility=full color=20,150,20
altColor=150,20,20 windowingFunction=mean
variableStep chrom=chr1
```

```
3003660 1
```

```
3004613 0.87804
```

```
3004636 0.86956
```

\*\*\*\*\*

**Bismark** and **BSMAP** is two popular BS-Seq align tools. This module can implement normalization for ".cov" files from Bismark(after v0.10.0) and ".bed" files from BSMAP.

Bismark: Normalization -m\_u 5,6 -i H1.cov -o out

BSMAP: Normalization --head -r\_t 4,5 -i H1.bed -o out

## (2) Identification of Methylation Patterns of regions

This module help you to get methylation patterns from the methylation files after data normalization.

-----  
USAGE: CpG\_MPs [options] [-i <folder>] [-o <folder>]  
-----

#### ARGUMENTS:

-h,--help,--man

Display the helpfile.

#### OPTION:

--methylation\_level\_U

The threshold of methylation level for unmethylated hotspot region.

The default methylation\_level\_U is set as 0.3.

--methylation\_level\_M

The threshold of methylation level for methylated hotspot region.

The default methylation\_level\_M is set as 0.7.

--hotspot\_CpG\_number

The successive CpG number (with methylation level  $\leq$  methylation\_level\_U or  $\geq$  methylation\_level\_M) threshold is used to find hotspot region (including unmethylated and methylated hotspot region).

The default hotspot\_CpG\_number is set as 3.

--abnormal\_CpG\_number

The CpG number with methylation level more than 0.5 is allowed to be  $\leq$  abnormal\_CpG\_number when extending the unmethylated hotspot region.

The CpG number with methylation level less than 0.5 is allowed to

be  $\leq$  abnormal\_CpG\_number when extending the methylated hotspot region. The default is set as 1.

```
--ucsc
    export extra ".bed" files in outdir to display on UCSC.

I/O:
-i,--indir
    The input methy folder (comes from
CpG_MPs_normalization).REQUIRED
-o,--outdir
    The output folder.REQUIRED
```

#### EXAMPLE:

```
CpG_MPs --ucsc -i NP_methy -o NP_pattern
    Get methylation patterns from NP_methy folder.The result is in NP_pattern.
    Export extra ".bed" files.
```

#### OUTPUT:

output will like this (chr13.patterns):

Chromosome	Start_position	End_position	Average_methylation_level	Variance	Methylation_pattern	Length
chr13	3000397	3127036	1	126640	0.897 0.120	
chr13	3127725	3357222	1	229498	0.913 0.107	
chr13	3357289	3360770	-1	3482	0.047 0.078	
chr13	3361137	3405393	1	44257	0.902 0.114	

The ucsc “.bed” format will like this:

```
track itemRgb="On" name="Methylation Region" description="methylation patterns regions"
visibility=3
```

chr13	3000436	3038679	MR	0	+	3000436 3038679 58,181,74
chr13	3038888	3127036	MR	0	+	3038888 3127036 58,181,74
chr13	3127725	3357354	MR	0	+	3127725 3357354 58,181,74
chr13	3357380	3360770	UMR	0	+	3357380 3360770 238,29,35
chr13	3360809	3405393	MR	0	+	3360809 3405393 58,181,74
chr13	3406351	3477301	MR	0	+	3406351 3477301 58,181,74
chr13	3477396	3477886	UMR	0	+	3477396 3477886 238,29,35
chr13	3478075	3536684	MR	0	+	3478075 3536684 58,181,74
chr13	3536764	3538770	UMR	0	+	3536764 3538770 238,29,35

### (3) Identification of Differentially and Conserved Methylated Regions

You can use this tool to get differentially methylation regions from mutiple methylomes.

```
USAGE: CpG_MPs_Dif_Con [options] [--patterns <folder_string>] [--methy
<folder_string>] [-o <folder>]
```

-----  
**ARGUMENTS:**

-h,--help,--man

Display the helpfile.

**OPTION:**

--V\_value

The value is set to estimate the identified differentially methylation region or conserved region covered by the extent of samples.

The identified region is more reliable with larger V\_value. The default V\_value is more than 0.5.

--CpG\_number

The identified differentially methylation region or conserved region must contain  $\geq$  CpG\_number CpG sites. The default CpG\_number is set as 4.

--ucsc

export extra ".bed" files in outdir to display on UCSC.

**I/O:**

--patterns

The patterns folder string.Split by ",".REQUIRED. like this:

NP\_patterns,ESC\_patterns

--methy

The input methy folder (comes from CpG\_MPs\_normalization).REQUIRED

Split by ",".like this NP\_methy,ESC\_methy

-o,--outdir

The output folder.REQUIRED

-----  
**EXAMPLE:**

CpG\_MPs\_Dif\_Con --ucsc --patterns NP\_patterns,ESC\_patterns,PGC\_patterns  
--methy NP\_methy,ESC\_methy,PGC\_methy -o NP\_ESC\_PGC

Get differentially methylation regions and conserved methylation regions from NP\_methy folder.

The result is in NP\_ESC\_PGC.Export extra ".bed" files.

**OUTPUT:**

output will like this (chr1.dmr):

Chromosome	Start_position	End_position	Length	Average_methylation_level	V_value
chr1	3943631	3943805	174	0.673_0.094_0.673	1
0.192695486142842					
chr1	4588643	4589123	480	0.233_0.706_0.233	1
0.308232831955781					
chr1	4679290	4679477	187	0.897_0.273_0.897	1
0.156921500047012					

chr1	4759803	4760445	642	0.561_0.261_0.561	1
0.694229364256945					
chr1	4774603	4774873	270	0.917_0.152_0.917	1
0.0759645439391204					

#### (4) Calculation the sequence features of the regions of DNA methylation patterns

This tool help you to get sequence feature for methylation patterns.

-----  
 USAGE: Seq\_Feature [--ref <folder>] [-i <folder>] [-o <folder>]  
 -----

#### ARGUMENTS:

--help,--man,-h  
                     display the helpfile.

#### I/O:

--ref  
                     the reference genome folder(.fa).REQUIRED  
 -i,--indir  
                     the input directory, methylation patterns folder. REQUIRED  
 -o,--outdir  
                     the output directory.REQUIRED  
 -----

#### EXAMPLE:

BRIEF(only required arguments):

Seq\_Feature --ref ref\_genome\_folder -i patterns/ -o seq\_feature/  
             get sequence feature for methylation patterns in patterns/.

#### OUTPUT:

the output will like this(chr1.seqf):

Chromosome	Start_position	End_position	Methylation_pattern	Length
Average_methylation_level	Variance	GC_content	CpG_OE	
chr1	3000574	3329885	1	329312
0.146981950106983			0.897	0.118
0.389521183558449			0.389521183558449	
chr1	3334584	3470553	1	135970
0.120034873515113			0.889	0.110
0.380260351548136			0.380260351548136	
chr1	3473945	3539119	1	65175
0.138771040623581			0.890	0.115
0.388937476026084			0.388937476026084	
chr1	3540334	3638079	1	97746
0.125952462499578			0.909	0.093
0.384077097784053			0.384077097784053	